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
- 1 Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. USDHHS. Centers for Disease Control, Atlanta.
- 2 Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, PA.
- 3 Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of health and humanservices, Centers for Disease Control and Prevention. Guideline for isolation precautions hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
- 4 U.S. Department of Health and Human Services. 1999. Biosafety in Microbiological and biomedical laboratories, HHS Publication (CDC), 4th.ed. U.S. Government Printing Office, Washington, D.C.
- 5 Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the Protection of Workers from Risks Related to Exposure to Biological Agents at Work (Seventh Individual Directive within The Meaning of Article 16 (1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

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